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GOAL OF THE EXPERIENCE

- 1) Measurement of the CCD of the DUT (poly-crystalline CVD diamond)
- 2) Comparison of the spectrum of a mip with a spectrum of highly ionising charge charged radiation

- Distribution of $Q(\text{measured})$ for mip and non-mip charged particles

- CCD Definition:
 - Charge Collection Distance:
 - $\text{CCD} = \text{track_length} \times Q(\text{measured})/Q(\text{generated})$
 - $\text{track_length} = \text{DUT thickness}$ for normally incident ionization radiation
 - $Q(\text{generated}) = |e| \times N(\text{ionizations})$
 - $Q(\text{measured}) = \text{induced signal calibrated in charge}$

SYNTHETIC DIAMOND (SINGLE CRYSTAL)

Property	Diamond	Silicon
Atomic number	6	14
Density, [g·cm ⁻³]	3.5	2.32
Band gap, [eV]	5.5	1.1
Resistivity, [Ω·cm]	> 10 ¹²	10 ⁵
Breakdown field, [V·μm ⁻¹]	1000 [107]	30 [108]
Electron mobility, [cm ² ·V ⁻¹ ·s ⁻¹]	1300~4500	1500
Hole mobility, [cm ² ·V ⁻¹ ·s ⁻¹]	2050~3800	500
Electron saturation velocity, [μm·ns ⁻¹]	141	100
Hole saturation velocity, [μm·ns ⁻¹]	96	100
Dielectric constant	5.6	11.7
Energy per e-h pair, [eV]	13~16 [111]	3.6
Av.min.ionizing signal per 100 m	3600	8000

Table 6.1: General properties of sCVD diamond compared with silicon in normal conditions

for a mip, average n. of ion pairs
produce in 100 micron

When the particles go through the diamond, the electron-hole pairs are generated due to the ionization processes. If a bias voltage is applied to the electrodes on the surfaces of the diamond, then an electric field is formed in the bulk of diamond. Under the force of this electric field, the electrons and holes start to move and this movement induces the current on the electrodes as described by the Shockley-Ramo theorem. Once the charge-carriers have arrived at the electrodes, the signals are “over” [115].


Figure 6.3: Charge generation and collection scheme

The charge collection distance (d_c) is considered as a common figure of merit to characterise the manufacture quality of a diamond, especially for the pCVD diamond, the d_c of which is usually much less than the thickness of the diamond (d). However, for the sCVD diamond, the d_c can sometimes be equal to or even greater than the thickness of the diamond.

As the average minimum-ionizing particle generates about 3600 electron-hole pairs in each 100 μm of diamond, the total charge generated (Q_{gen}) can be calculated as $Q_{gen} = 36 \times d$, where d is in μm scale.

If the electrons and holes get trapped or recombined after separating by the distance d_c , then the charge collected by the electrodes (Q_{coll}) can be expressed as:

$$Q_{coll} = Q_{gen} \times \frac{d_c}{d} \quad (6.2)$$

Knowing the charge collected, the charge collection distance can be calculated as:

for a mip, average n. of ion pairs
produce in 1 micron

$$d_c = \frac{Q_{coll}}{36 \text{ e}^- / \mu\text{m}} \quad (6.3)$$

The ratio between the total generated charge and the charge collected by the electrodes is defined as the charge collection efficiency (CCE):

$$CCE = \frac{Q_{coll}}{Q_{gen}} \quad (6.4)$$


Combining Eq. 6.2 and Eq. 6.4, d_c can then be written as:

$$d_c = d \times CCE \quad (6.5)$$

The CCE has a dependence on the applied bias voltage and it can be measured using the beta sources as described in Section 6.2.1.

HOW CONCEPTUALLY


- let's assume to have mono-energetic ionising particles ($E > E_{\text{mpv}}$) traversing the ideal diamond normally
- one can estimate, using data from literature, the $\langle \Delta E \rangle$ lost in a 500 μm path in diamond (3600 ionizations/100 μm)
- the measured E_{loss} will fluctuate with a Landau-Vavilov distribution: m.p. E_{loss} (mpv = [1]) and FWHM (related to the second parameter [2]) will parameterize the shape of the distribution (in TMath::Landau)
- $\langle \Delta E \rangle$ or mpv can be related to the average or mpv of the signal charge integral or signal pulse height distributions (amplifier gain and noise and readout input impedance define how the energy loss and the signal parameter are related)
- the signal parameter distribution can be measured on the DUT with the same readout chain (same calibration) => the known E_{loss} in the ideal detector, provides the absolute calibration.

HOW TO EMULATE A M.I.P. ?

- One step back:
 - how to select a sample of ionizing particles releasing (in average) always the same amount of energy ?
 - no need of monoenergetic particles, a mip is enough
 - for $E > E_{\text{mip}}$ average Eloss is very stable
- muons from cosmic rays ?
 - too low rate
 - a not too well defined sample (E distribution)


Energy loss of electrons in C

- $E(\text{kinetic}) > 800\text{KeV} \Rightarrow e^-$ is a m.i.p. in Carbon/Diamond/Graphyte

ELECTRON ENERGY LOSS IN CARBON

- zoom on the region **800 KeV - 2 MeV**
 - no appreciable variation of the mean E_{Loss}

HOW TO EMULATE A M.I.P. ?

<https://www.nist.gov/pml/productsservices/physical-reference-data>

- radiation/dosimetry data

<https://www.nist.gov/pml/stopping-power-range-tables-electrons-protons-and-helium-ions>

Contents:


1. [Introduction](#)
2. [ESTAR](#): Stopping Powers and Ranges for Electrons
3. [PSTAR and ASTAR](#): for Protons and Helium Ions (alpha particles)
4. [References](#)
5. [Appendix](#): Significance of Calculated Quantities


Access the Data:

1. [Electrons](#)
2. [Protons](#)
3. [Helium Ions](#)

- $E(\text{kinetic}) > 800\text{KeV} \Rightarrow e^-$ is a m.i.p. in Carbon/Diamond/Graphyte

Energy loss of electrons in C

The ESTAR program calculates stopping power, density effect parameters, range, and radiation yield tables for electrons in various materials. Select a material and enter the desired energies or use the default energies. Energies are specified in MeV, and must be in the range from 0.001 MeV to 10000 MeV.

[Help](#)

[Text version](#)

[Material composition data](#)

Select a common material:

6: Carbon, Amorphous (density 2.0g/cm³)

or enter a [unique material](#)

- Graph stopping power:**
 - Total Stopping Power
 - Collision Stopping Power
 - Radiative Stopping Power

Graph density effect parameter

Graph CSDA range

Graph radiation yield

No graph

Additional Energies (optional):

Use energies from a file*

Choose File no file selected

or

Use energies entered below (one per line)

Include default energies

Note: Only stopping powers and the density effect parameter will be calculated if additional energies are used.

Submit

Reset


ELECTRON SOURCES: RADIOACTIVE SOURCE

- ^{90}Sr is a product of nuclear fission. It is present in significant amount in spent nuclear fuel and in radioactive waste from nuclear reactors and in nuclear fallout from nuclear tests. For thermal neutron fission as in today's nuclear power plants, the fission product yield from U-235 is 5.7%, from U-233 6.6%, but from Pu-239 only 2.0%.
- ^{90}Sr undergoes β^- decay with a half-life of 28.79 years and a decay energy of 0.546 MeV distributed to an electron, an anti-neutrino, and the yttrium isotope ^{90}Y , which in turn undergoes β^- decay with half-life of 64 hours and decay energy 2.28 MeV distributed to an electron, an anti-neutrino, and ^{90}Zr (zirconium), which is stable. *Note that $^{90}\text{Sr}/\text{Y}$ is almost a pure beta particle source; the gamma photon emission from the decay of ^{90}Y is so infrequent that it can normally be ignored.*

Strontium-90	
Full table	
General	
Name, symbol	Strontium-90, ^{90}Sr
Neutrons	52
Protons	38
Nuclide data	
Natural abundance	syn
Half-life	28.79 years
Decay products	^{90}Y
Decay mode	Decay energy
Beta decay	0.546 MeV

SPECTRUM OF ELECTRONS FROM BETA DECAYS


- max E_{electron} given by the Q value of the decay
- neutrino E is $Q - E_{\text{electron}}$

Beta spectrum : a sharing of energy:

In a beta decay like that of a bismuth-210 nucleus, the decay energy is shared between the nucleus, an electron and an antineutrino. As the nucleus, is much more massive than the two others (its mass is 320 000 times the mass of the electron) it takes away a negligible amount of kinetic energy. The energy is therefore shared between the electron and the antineutrino. The antineutrino escaping detection, one observe only the beta electron with a variable energy. The figure display the characteristic energy repartition of beta electrons – the beta spectrum – of the bismuth-decay

©IN2P3

Examples of average beta energies

Since the energy of beta electrons emitted by a radionuclide. is not unique, one compares generally radionuclides by their average beta energies. These energies vary considerably : for instance, the average energy of beta emitted by tritium is a hundred times lower than that of phosphorus-32. The average energies are below 1 MeV, much lower than those of alpha particles (usually above 4 MeV).

©IN2P3

The prevalence of low energy beta electrons is beneficial for radiological protection, because low energy electrons are easier to stop. In radiological protection, one is more interested in the average energy of electrons than by their maximum energy. This average energy varies within wide proportions, ranging from 5.69 keV in the case of tritium to 695 keV in the case of phosphorus-32 a powerful beta emitter.

Compared to alpha particles, the energies of beta electrons are much lower, below 1 MeV in most of the cases while the energies of alpha particles are always above 4000 keV (4 MeV). The lifetimes (half-lives) are much shorter, with the exception of potassium-40.

A beta decay often comes with the emission of a gamma-ray due to the desexcitation of the nucleus. This emission decreases the energy to be shared between the electron and the antineutrino. For instance, the available energy in the caesium-137 beta decay is 1176 keV, but in 95 % of the cases it is accompanied by the emission of a characteristic 662 keV gamma in which case the available energy is only 514 keV. The observed spectrum beta becomes the sum, in the proportions of 5% and 95%, of the two spectra corresponding to the modes without or with gamma.

ELECTRONS FROM A ^{90}Sr SOURCE

- The decay energy Q is the mass difference Δm between the parent and the daughter atom and particles
- beta- decay:
 - $n \rightarrow p e^- \text{ anti-}\nu_e$
 - $Q = 939.565 - 938.272 - 0.511 = 0.782 \text{ MeV}$
 - $^{90}\text{Sr} \rightarrow ^{90}\text{Y} e^- \text{ anti-}\nu_e$
 - $Q=546 \text{ MeV} \rightarrow$ the maximum energy of electrons from a ^{90}Sr decay is 0.546 MeV smaller than E_{mip} ;
 - the $\langle E \rangle$ of e^- is 196 KeV
 - $^{90}\text{Y} \rightarrow ^{90}\text{Zr} e^- \text{ anti-}\nu_e$
 - $Q = 2.28 \text{ MeV}$
 - the $\langle E \rangle$ of e^- is 1 MeV

Main isotopes of yttrium

iso	NA	half-life	DM	DE (MeV)	DP
^{87}Y	syn	3.4 d	ϵ	–	^{87}Sr
			γ	0.48, 0.38D	–
^{88}Y	syn	106.6 d	ϵ	–	^{88}Sr
			γ	1.83, 0.89	–
^{89}Y	100%	is stable with 50 neutrons			
^{90}Y	syn <i>Z=39, N=51</i>	2.7 d	β^-	2.28	^{90}Zr
			γ	2.18	–
^{91}Y	syn	58.5 d	β^-	1.54	^{91}Zr
			γ	1.20	–


Standard atomic weight (A_r) 88.905 84(2)^[1]

[view](#) · [talk](#) · [edit](#)

Zirconium
stable
Z=40, N=50


Zirc

ELECTRON SPECTRUM IN ^{90}Sr

- after cutting electrons with energies < 0.8 MeV in the e- spectrum of a ^{90}Sr source, we are left with a m.i.p. sample

ELECTRON SPECTRUM IN ^{90}Sr

- $E(\text{kinetic}) > 800\text{KeV} \Rightarrow e^-$ is a m.i.p. in Carbon/Diamond/Graphyte



Figure 6.2: Left: the mean energy loss of an electron in diamond; Right: Landau distribution simulated using Geant4

- https://groups.lal.in2p3.fr/atf2/files/2015/05/Chapter6_v3.pdf

SETUP: MEASUREMENT MODE


Shielding case

Radioactive source

Aluminum plate

DUT pCVD diamond
calibration sCVD diamond
energy filter
collimator
triggering diamond


SETUP: CALIBRATION MODE



DIAMOND

- from ESTAR db

NIST_stoppingPower_e_amorfousCarbon.dat

NIST_stoppingPower_e_amorfousCarbon.dat

BRASS

- Brass (60% Cu 40% Zn);
- density = 8.5 g/cm³;
- <http://physics.nist.gov/PhysRefData/St>

Fill out the following form to define your unique material, or return to the [common material form](#).

[Help](#)

[Text version](#)

Material Name:



Density: (g/cm³)

**Enter the formulae and relative weights separated by a space for each item.
One compound or element per line. For example, salt water with 10 % (by weight) salt:**

H2O 0.9
NaCl 0.1

Note: Weights not summing to 1 will be normalized.

NIST_stoppingPower_e_Brass_60_40_8500.dat

CAL. DIAMOND & DUT WITH BRASS 0.05CM

Diamond

Density:

3.5 (g/cm³)

Mean Excitation Energy:

81.0 (eV)


Atomic Constituents:

C

1.

Stopping power in Diamond (3.5g/cm³) and Brass60-40 (8.5g/cm³)

StoppingPowerIntegrateIdEdx.C

Estart 2.3 MeV ... very good m.i.p.


Elost in DUT 0.27925 MeV

in Calib Diamond 0.277849 MeV

in Brass 0.560274 MeV Brass thickness = 0.05 cm

in Triggering Diamond 0.278285 MeV

Energy loss/micron vs depth [KeV]


electron survives the materials of the setup; final energy = 0.90 MeV

1MeV e are absorbed in brass, no trigger


NO BRASS, CAL. DIAMOND & DUT

Stopping power in Diamond (3.5g/cm³) and Brass60-40 (8.5g/cm³)


Estart 1 MeV
 Elost in DUT 0.27925 MeV
 in Calib Diamond 0.277849 MeV
 in Brass 0. MeV Brass thickness = 0 cm
 in Triggering Diamond 0.459445 MeV
 Effective track length in triggering Diamond = 0.044 cm
 Estart 1 MeV
 Estart - ElostinDUT 0.721 MeV
 Estart - ElostinDUT - ElostinCD **0.443 MeV**



Estart 2.3 MeV
 Elost in DUT 0.27925 MeV
 in Calib Diamond 0 MeV
 in Brass 0.277849 MeV Brass thickness = 0 cm
 in Triggering Diamond 0.275586 MeV
 Energy after triggering Diamond = 1.46732 MeV

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